

ESPRIA and Media: From a hearing system to a multimedia system.

Connecting to today's world is very simple: ESPRIA and Media bring you excellent sound quality and optimum levels of comfort in all modern media – using Bluetooth!

- Connect to your cell/phone
- Connect to your television (with Media Transmitter)
- Connect to external signal sources, e.g. your MP3 player or CD and DVD player from your home entertainment system.

Classically beautiful and refined: eMote2.

ESPRIA can, of course, also be used with a classic remote control.




At a glance: ESPRIA by HANSATON.

- ESPRIA is the hearing system for a new generation that opens up more possibilities and perspectives to modern people.
- ESPRIA is the innovative HANSATON range of hearing systems that combines the latest technology with excellent design in five models of your choice.
- ESPRIA is an excellent choice; from the individual ITE hearing system through the award-winning X-Mini and tiny Mini hearing systems to the elegant ultimate performances in the Slim and Power format.
- ESPRIA is the new Business Class manufactured by HANSATON – premium technology that cannot fail to convince:
 - i-com2 continuously and automatically synchronizes both hearing systems for improved speech intelligibility and directional hearing, especially in a noisy environment.
 - Situation Optimizer, the all-rounder that adjusts immediately to every situation.
 - Active Feedback Block with the latest phase cancellation feedback management.
 - Media, the elegant and innovative remote control for top sound quality and optimum comfort – linked via Bluetooth to all modern media.

ESPRIA is the perfect hearing system for a new generation of modern people.

Your hearing instrument professional will be happy to help you:

Enjoying the freedom of living life to the full.



"We encounter new perspectives,

more possibilities and innovative technologies every day. New ideas that simplify and enrich our lives are continually being conceived. And we accept them with curiosity and enthusiasm. The latest cell phone, the ideal home entertainment system, the most environmentally-friendly car. And now: the optimal hearing system. ESPRIA. Hearing. Communicating. Being part of it. My new freedom."



ESPRIA by HANSATON:
Five different versions for your
individual desires and needs.



ESPRIA In-The-Ear hearing systems:
valuable and unique natural designs.

ESPRIA X-Mini – The award-winning designer housing is so tiny it is almost invisible.



The new **ESPRIA Mini**, Business Class technology with an elegantly-styled look.

ESPRIA Slim. Exceptional:
high performance with
extra-slim design.



NEW – available
from Nov. 09




ESPRIA Power, the system that will
meet even the most rigorous standards.

"Highly individual: Nature and technology working together."

HANSATON the hearing system manufactory uses the individual mold of your auditory canal to create your personalized In-The-Ear profile – with digital precision in just a few days. Whether it is the CIC, Mini-Canal, Canal or Concha model – all ESPRIA In-The-Ear hearing systems fit comfortably, securely and precisely in your ear.




CIC



Mini-Canal




Concha



Canal

Welcome to the latest generation of ITE hearing systems.

- Maximum venting (air channel system for optimum ventilation) for unrivaled wearing comfort.
- Anti-cerumen concepts to protect the sensitive technology.
- Express digital manufacturing service.
- Suitable for mild to severe hearing losses.
- With right/left side identification.



Little works of art: the individually handmade „Legance“ surface emulates your natural skin structure so that it is hard to tell the difference.

"Two technologies, one effect: Minimalism from HANSATON."

A tiny housing holds a completely new attitude to life. ESPRIA X-Mini and ESPRIA Mini from HANSATON are among the smallest Behind-The-Ear hearing systems in the world.

Fascinating technology: experience the ESPRIA X-Mini.

- The latest technology and removal of the receiver ensure that the housing is tiny.
- Modern openness: natural sound quality without blocked ears.
- Three receivers for all hearing needs (45/55/65 dB) and in several lengths.
- Suitable for mild to profound hearing losses.
- Fitting with a range of domes or earmolds.
- Right/left side identification and on/off with battery compartment.

Perfect design in seven natural color tones.



The sophisticated 2-tone special paint makes them look even slimmer and the replaceable shell concept enhances their attractiveness: your hearing instrument professional can simply fit a new top shell.

A tiny design piece: discover the ESPRIA Mini.

- Suitable for mild to moderate hearing losses.
- Fitting with a mini-tube or earmold.
- No operating controls on the hearing system.
- Nano coating to repel dirt and moisture.
- Right/left side identification and on/off with battery compartment.



Awarded with



and nominated for the design prize
of the Federal Republic of Germany.



"Range also includes strong partners for more intensive requirements."

ESPRIA relies on the Slim and Power models if more support is required. The most modern Behind-The-Ear hearing systems with an optical design and sophisticated color palette that conceals the actual amount of power that they offer.

ESPRIA Slim: the ultimate elegant performance.

- Suitable for mild to moderate hearing losses.
- Fitting with an earmold (mini-tube is possible).
- Patented push scroll – volume and program are set with an operating control.
- Replaceable shell concept and 2-tone special paint.
- Right/left side identification and on/off with battery compartment.



ESPRIA Power: beautiful yet strong.

- Up to 70 dB gain with medium to profound hearing losses.
- Can be fitted with an earmold.
- Simple operating controls for volume and program.
- Battery safety lock, audio adapter and accessories.
- Right/left side identification and on/off with battery compartment.



Also in five cheerful children's colors.



The latest technology with intelligent management.

A complete all-rounder: the Situation Optimizer.

Listening to music, holding a conversation, concentrating on just one sound source even in noisy surroundings: ESPRIA continuously and automatically analyzes your environment and active requirements. Enjoy hearing that is so perfectly managed that you will no longer even notice your hearing system.

Only let the good stuff through: the Active Feedback Block.


ESPRIA uses the latest generation of phase cancellation feedback management. A highly-technical system ensures that feedback is recognized immediately and that not a single sound is allowed to reach you.

Intelligent communication: two in one with i-com2*!

Two hearing instruments that communicate with each other and automatically adjust to the other: i-com2 provides a completely synchronized hearing experience with balanced sound and excellent speech intelligibility – even in a noisy environment.



eMote2



Media 

*A non-wireless version of ESPRIA is available in some countries, only.
The i-com2 function is not available in these countries.