

## At a glance: LOONA by HANSATON.

- LOONA offers exclusive design and the latest technology combined with attractive conditions from HANSATON.
- LOONA is an exceptional hearing system range in the Comfort Class with five models of your choice to meet a range of personal preferences and any possible requirements.
- LOONA is an excellent choice; from the individual ITE hearing system through the award-winning X-Mini and tiny Mini hearing systems to the elegant ultimate performances in the Slim and Power format.
- LOONA comprises intelligent technology for modern customers who value LOONA for the following highlights, to name just a few:
  - i-com2 continuously and automatically synchronizes both hearing systems for extremely improved speech intelligibility and directional hearing, especially in a noisy environment.
  - Active Feedback Block with the latest phase cancellation feedback management.
  - Sound Impulse Management for more pleasant hearing without interruptions from sudden pulse-like signals.

LOONA is the hearing system that makes you feel secure in knowing you have made the right choice in terms of design, price and performance.

Your hearing instrument professional will be happy to help you:



Hansaton Akustik GmbH  
 Stueckenstrasse 48 · 22081 Hamburg · Germany  
[www.hansaton.com](http://www.hansaton.com)



Definitely  
the right  
decision.

loona 

The unbeaten combination  
of price and performance.



# "It's clear to us very quickly – say yes to LOONA!"

HANSATON makes it easy to make the right choice, because the first-class quality provided by the manufacturer of modern hearing systems is available for every desire and requirement. In LOONA, HANSATON has combined the highest levels of attractiveness from every area: design, performance and price. As if it were made for us."



Get to know the LOONA family: the optimum hearing experience in five different versions of your choice.



**LOONA In-The-Ear hearing systems:** individually fitted and almost invisible.

**LOONA X-Mini** – tiny multi award winning design housing with an external receiver.



The new **LOONA Mini:** the outside is lovingly designed in a refined style while the inside offers the latest technology.

**LOONA Slim,** award-winning design with extra slim looks.



NEW – available from Nov. 09



**LOONA Power,** the hearing system for particularly intensive requirements.

# “Made-to-measure – hearing systems designed just for me.”

HANSATON the hearing system manufactory uses the individual mold of your auditory canal to create your personalized In-The-Ear profile – with digital precision in just a few days. Whether it is the CIC, Mini-Canal, Canal or Concha model – all LOONA In-The-Ear hearing systems fit comfortably, securely and precisely in your ear.



CIC



Mini-Canal




Concha



Canal

## Welcome to the latest generation of ITE hearing systems.

- Maximum venting (air channel system for optimum ventilation) for unrivaled wearing comfort.
- Anti-cerumen concepts to protect the sensitive technology.
- Express digital manufacturing service.
- Suitable for mild to severe hearing losses.
- With right/left side identification.



Little works of art: the individually handmade „Legance“ surface emulates your natural skin structure so that it is hard to tell the difference.

# "Perfect technology in a tiny designer casing."

A tiny housing holds a completely new attitude to life. LOONA X-Mini and LOONA Mini from HANSATON are among the smallest Behind-The-Ear hearing systems in the world.

## Fascinating technology: experience the LOONA X-Mini.

- The latest technology and removal of the receiver ensure that the housing is tiny.
- Modern openness: natural sound quality without blocked ears.
- Three receivers for all hearing needs (45/55/65 dB) and in several lengths.
- Suitable for mild to profound hearing losses.
- Fitting with a range of domes or earmolds.
- Right/left side identification and on/off with battery compartment.

## Perfect design in seven natural color tones.



The sophisticated 2-tone special paint makes them look even slimmer and the replaceable shell concept enhances their attractiveness: your hearing instrument professional can simply fit a new top shell.

## A tiny design piece: discover the LOONA Mini.

- Suitable for mild to moderate hearing losses.
- Fitting with a mini-tube or earmold.
- No operating controls on the hearing system.
- Nano coating to repel dirt and moisture.
- Right/left side identification and on/off with battery compartment.



Awarded with



and nominated for the design prize  
of the Federal Republic of Germany.



# "For more intensive requirements – solutions with more power."

LOONA relies on the Slim and Power models if more support is required. The most modern Behind-The-Ear hearing systems with an optical design and sophisticated color palette that conceals the actual amount of power that they offer.

## LOONA Slim: the ultimate elegant performance.

- Suitable for mild to moderate hearing losses.
- Fitting with an earmold (mini-tube is possible).
- Patented push scroll – volume and program are set with an operating control.
- Replaceable shell concept and 2-tone special paint.
- Right/left side identification and on/off with battery compartment.



## LOONA Power: beautiful yet strong.

- Up to 70 dB gain with medium to profound hearing losses.
- Can be fitted with an earmold.
- Simple operating controls for volume and program.
- Battery safety lock, audio adapter and accessories.
- Right/left side identification and on/off with battery compartment.



## Also in five cheerful children's colors.



# The latest technology that cannot fail to convince.

## **Intelligent communication: two in one with i-com2\*!**

Two hearing instruments that communicate with each other and automatically adjust to the other: i-com2 provides a completely synchronized hearing experience with balanced sound and excellent speech intelligibility – even in a noisy environment.

## **Only let the good stuff through: the Active Feedback Block.**

LOONA uses the latest generation of phase cancellation feedback management. A highly-technical system ensures that feedback is recognized immediately and that not a single sound is allowed to reach you.

## **Improved hearing with no disturbances: thanks to Sound Impulse Management.**

Pulse-like signals e.g. dishes rattling, paper rustling and doors slamming are recognized and reduced effectively. Of course, this is done without reducing hearing quality at the same time – just intelligent sound management!

## **Classically beautiful and refined: eMote2\*.**

LOONA can, of course, also be used with a remote control. Simple to use with a discreet design – eMote2 includes the following functions: Adjusting the volume, switching hearing programs, switching the hearing systems on/off.

